

Reading a Triple Beam Balance

<http://www.ohaus.com/products/education/weblab/TBBread.html>

1. To begin, click "Next".
2. Turn the sound off on your laptop.
3. Click on "New Unknown".
4. Read the TBB to the nearest 10th of a gram.
5. Click on "Check answer".
6. Record your data in the data table.



Cut along dotted lines and paste into lab journal

Riders					
#	Hundreds	Tens	Ones	Reading	Actual
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					



Conclusion: 2-3 sentences on what you learned